
*Scan QR Code with your WebCam for downloading the user instruction manual.*

# PR Series Balances

## Quick Start Guide

EN

ES

FR

*Please download the user manual from [www.ohaus.com](http://www.ohaus.com).*

## 1. SAFETY INFORMATION

This manual contains installation, operation and maintenance instructions for PR Series balance. Please read this manual completely before installation and operation.

### Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

<b>WARNING</b>	For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.
<b>CAUTION</b>	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided.
<b>Attention Note</b>	For important information about the product For useful information about the product

### Warning Symbols

General Hazard

Electrical Shock Hazard

Alternating Current

Direct Current

### Safety Precautions

**CAUTION:** Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Verify that the AC adapter's input voltage range and plug type are compatible with the local AC main power supply.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Do not position the balance such that it is difficult to reach the power connection.
- The balance is for indoor use only. Do not operate the equipment in hazardous or unstable environments.
- Operate the equipment only under ambient conditions specified in these instructions.
- Do not drop loads on the pan.
- Use the balance only in dry locations.
- Disconnect the equipment from the power supply when cleaning.
- Use only approved accessories and peripherals.
- Service should only be performed by authorized personnel.

### Intended Use

Use the instrument exclusively for weighing as described in the operating instructions. Any other type of use and operation beyond the limits of technical specifications without written consent from OHAUS, is considered as not intended.



This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use.

If the instrument is not used according to these operating instructions, the intended protection of the instrument may be compromised and OHAUS assumes no liability.

## 2. INSTALLATION

### 2.1 Select the Location

Avoid heat sources, rapid temperature changes, air current or excessive vibrations. Allow sufficient space.


### 2.2 Leveling

Be sure the balance is level before it is used or after its location is changed.

The PR balance has a level bubble in a small round window beside the display.

To level the balance, adjust the 2 leveling feet until the bubble is centered in the circle.

Please refer to the right figure for leveling.

### 2.3 Connecting Power and Acclimatising the Balance

Connect the DC output connector to the power receptacle on the rear of the balance. Then connect the AC adapter plug to a suitable electrical outlet.

#### Acclimatising

It is suggested that the balance should not be used until it has been connected to power and acclimatised to the environment for a certain period of time. In the case of a balance with the precision above 0.1 mg, the acclimatisation time should be 1.5 hours; in the case of balance with the precision of 0.01 mg, the acclimatisation time should be more than 4 hours.

### 2.4 Connecting the Interface

The PR balance has a RS232 port.

Use the RS-232 port to connect either to a computer or a printer with a standard (straight-through) serial cable.

Interface connections on the rear of the balance

RS232

RS232 : Used to connect to PC or Printer

Note: See the Printing section for Connecting, Configuring and Testing the Printer / Computer Interface.


## 2.5 Initial Calibration

For best results, the instrument's calibration should be checked prior to first use. If adjustment is needed, refer to the Calibration section of the instruction manual.

## 3. OPERATION

### 3.1 Overview of Controls and Display

#### CONTROLS

#### CONTROL FUNCTIONS

Button	 Zero Yes	 Print Unit No	 Function Mode Back	 Tare Menu Cal Exit
<b>Primary Function (Short Press)</b> 	<b>On / Zero</b> <ul style="list-style-type: none"> <li>If the balance is Off, turns on the balance.</li> <li>If the balance is On, sets zero.</li> </ul>	<b>Print</b> <ul style="list-style-type: none"> <li>Sends the current displayed value to the serial interface.</li> </ul>	<b>Function</b> <ul style="list-style-type: none"> <li>Operation is dependent on the application mode.</li> </ul>	<b>Tare</b> <ul style="list-style-type: none"> <li>Performs tare operation.</li> </ul>
<b>Secondary Function (Press and Hold)</b> 	<b>Off</b> <ul style="list-style-type: none"> <li>Zeros the current value.</li> </ul>	<b>Unit</b> <ul style="list-style-type: none"> <li>Changes weighing units.</li> </ul>	<b>Mode</b> <ul style="list-style-type: none"> <li>Changes application mode.</li> </ul>	<b>Menu-Cal</b> <ul style="list-style-type: none"> <li>Enters the main menu. Calibration is the first sub-menu.</li> <li>Views the preset Tare value.</li> </ul>
<b>Menu Function (Short Press)</b> 	<b>Yes</b> <ul style="list-style-type: none"> <li>Accepts the current (blinking) setting on the display.</li> </ul>	<b>No</b> <ul style="list-style-type: none"> <li>Rejects the current (blinking) setting on the display.</li> <li>Increments a value being entered.</li> </ul>	<b>Back</b> <ul style="list-style-type: none"> <li>Reverts back to previous menu item.</li> <li>Decrements a value being entered.</li> </ul>	<b>Exit</b> <ul style="list-style-type: none"> <li>Immediately exits the sub-menu.</li> <li>Aborts a calibration in progress.</li> </ul>

DISPLAY

<p>Net (NET) Stability (*) Negative (-) Centre of zero (&gt;0&lt;) Pre-tare (PT)</p>		<p>Result Field: depending on application</p> <p>Unit</p>
--	--	---

3.2 Principal Functions

- Weighing:** First press Zero to set the display to zero. Place an object on the pan. The display indicates the gross weight.
- Taring:** With no load on the pan, press Zero to set the display to zero. Place an empty container on the pan and press Tare. Add objects to the container and its net weight is displayed. After the container and the objects are removed, the load will be displayed as a negative number. Press Tare to clear.
- Zero:** Press Zero to zero the balance.

4. MAINTENANCE

4.1 Cleaning

**WARNING:** Disconnect the balance from the power supply before cleaning. Make sure that no liquid enters the interior of the balance.

Clean the balance at regular intervals.  
Housing surfaces may be cleaned with a lint-free cloth slightly dampened with water or a mild cleaning agent.  
Glass surfaces may be cleaned with a commercial glass cleaner.  
**Attention:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

4.2 Troubleshooting

If the troubleshooting section does not resolve your problem, contact an Authorized Ohaus Service Agent. Please visit our website [www.ohaus.com](http://www.ohaus.com) to locate the Ohaus office nearest you.

## 5. TECHNICAL DATA

The technical data is valid under the following conditions:

- Indoor use only
- Powered by external approved LPS AC/DC adaptor
- Altitude: Up to 2000 m
- Operating temperature range: 5 °C to 40 °C
- Specified Temperature range: 10°C to 30°C
- Humidity: maximum relative humidity 80% for temperatures up to 30°C, decreasing linearly to 50% relative humidity at 40°C
- Mains supply voltage fluctuations: up to  $\pm 10\%$  of the nominal voltage
- Installation category II
- Pollution degree: 2

Model	PR124	PR224			
	PR124/E	PR224/E	PR223/E	PR423/E	PR523/E
Capacity	120	220	220	420	520
Readability d (g)	0.0001	0.0001	0.001	0.001	0.001
Repeatability (std. dev.) (g)	0.0001	0.0001	0.001	0.001	0.001
Linearity (g)	$\pm 0.0002$	$\pm 0.0002$	$\pm 0.002$	$\pm 0.002$	$\pm 0.002$
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A				

Model	PR1602/E	PR2202/E	PR4202/E	PR4201/E	PR6201/E
Capacity	1600	2200	4200	4200	6200
Readability d (g)	0.01	0.01	0.01	0.1	0.1
Repeatability (std. dev.) (g)	0.01	0.01	0.01	0.1	0.1
Linearity (g)	$\pm 0.02$	$\pm 0.02$	$\pm 0.02$	$\pm 0.2$	$\pm 0.2$
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A				

Model		PR523N				
	PR323N/E	PR523N/E	PR322N/E	PR522N/E	PR822N/E	PR2202N/E
Class	II	II	II	II	III	II
Capacity	320	520	320	520	820	2200
Verification scale interval e (g)	0.01	0.01	0.01	0.01	0.1	0.1
Readability d (g)	0.001	0.001	0.01	0.01	0.1	0.01
Repeatability (std. dev.) (g)	0.001	0.001	0.01	0.01	0.01	0.01
Linearity (g)	$\pm 0.002$	$\pm 0.002$	$\pm 0.02$	$\pm 0.02$	$\pm 0.02$	$\pm 0.02$
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A					

## EN-6

Model	PR4202N/E	PR5202N/E	PR2201N/E	PR4201N/E	PR5201N/E	PR6201N/E
Class	II	II	II	II	II	III
Capacity	4200	5200	2200	4200	5200	6200
Verification scale interval e (g)	0.1	0.1	0.1	0.1	0.1	1
Readability d (g)	0.01	0.01	0.1	0.1	0.1	1
Repeatability (std.dev.) (g)	0.01	0.01	0.1	0.1	0.1	0.1
Linearity (g)	±0.02	±0.02	±0.2	±0.2	±0.2	±0.2
Power Supply	Power input: 100-240V ~ 200mA 50-60Hz 12-18VA Power output: 12 VDC 0.5A					

## 6. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a> .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1 UL Std. No. 61010-1

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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P/N 30431476B © 2018 Ohaus Corporation, all rights reserved.
Ohaus Corporation