



MEASUREMENT & CALIBRATION TECHNOLOGIES DIVISION

8600 Somerset Drive, Largo, Florida 33773 U.S.A.

Telephone: 800-527-9999 Fax: 727-539-6882

Email: scott.blanchard@AMETEK.com

Web www.ametektest.com

1-1-2020

Subject:

Minimum PC requirements for running Nexygen Plus Version 4

Dear Customer,

PC Requirements for running our Lloyd Instrument Nexygen Plus software are below:

- PC with i3 N5000 processor or better @ 2.7 GHz
- 64 bit operating system
- 4GB of RAM or more
- 2 Available USB Ports
- Microsoft Office 2016. Including Access, Word, and xlsx for reporting and security features.
- Windows 10 Operating System

Best Regards,

Name: Scott Blanchard

Position: Sales Manager

AMETEK Measurement & Calibration Technologies Division

CHATILLON, DAVENPORT, LLOYD INSTRUMENTS Brands

Email: scott.blanchard@ametek.com

Phone: 518-605-7891

Web www.ametektest.com