



		GP Series	GX-L Series
Dimensions	Main unit	Swing arm type: 14.65" (W) x 24.2" (D) x 5.1" (H) Separate display type: 13.5" (W) x 17.5" (D) x 5.1" (H) Display: 10.25" (W) x 6.5" (D) x 7.9" (H)	
	Weighing pan	15.1" x 13.5"	
Weight	Main unit	35 lb	
Material	Main unit	Aluminum die-casting (Display upper case: ABS)	
	Weighing pan	SUS316	
	Display protective cover	PET	
Weighing	Counting	✓	✓
	Comparator	3 levels	3 levels/5 levels
	Animal weighing	✓	✓
	Statistical calculation	Addition only	✓
	Density measurement	✓	✓
	Flow rate detection (FRD)	x	✓
	Impact shock detection (ISD)	x	✓
Function	Internal weight sensitivity adjustment	✓	✓
	Automatic self sensitivity adjustment	✓	✓
	Buzzer	x	✓
	Repeatability test	✓	✓
	Automatic precision assessment	x	✓
	Minimum weight alert	x	✓
	Data memory Weight value	✓	✓
	Data memory Sensitivity adjustment log	✓	✓
	Data memory Unit weight value	✓	✓
	Data memory Comparator value	✓	✓
	Data memory Tare value	✓	✓
	Comparator enlarged display	✓	✓
	Password lock & User Login	x	✓
	Spirit level LED	x	✓
	Automatic level detection & warning	x	✓
Output	GLP output	✓	✓
	Date & Time output	✓	✓
	Interval output	x	✓
Interface	RS-232C	✓	✓
	USB (Quick USB)	x	✓
	USB (Virtual COM)	x	✓
	Comparator output	Optional (GP-04)	Optional (GXL-04)
	Current loop output	Optional (GP-04/GP-06)	x
	Analog output (voltage)	Optional (GP-06)	Optional (GXL-06)
	LAN	No	Optional (GXL-08)
	Bluetooth	No	Optional (GXL-27)
	External key input	RS-232C	Optional (GXL-04)
	Underhook	Optional (GP-21/22)	Optional (GP-21/22)
Others	Animal weighing pan	Optional (GP-12)	Optional (GP-12)
	Waterproof RS-232C cable (5m)	x	AX-KO2737-500